

arkray

Advanced Blood Glucose Test Meter

GLUCOCARD X-mini plus

GT-1960



More advanced SMBG meter
which successfully benefits
the user with no coding,
top specifications,
and stylish design.

arkray, inc.

Advanced Blood Glucose Test Meter

GLUCOCARD X-mini plus

GT-1960



Design

Impressive & Stylish design

Operation

Comprehensive display
New user friendly functions

Small Sample

0.3uL micro sampling

Maximum flexibility

Auto coding
Short fill detection
Strip discard button
Test strip port light
Backlight LCD
External output



Test strip for GLUCOCARD X-mini plus

Specification			
Test	Glucose in whole blood	Memory	500 test results
Sample Volume	Approx. 0.3uL	Operating Environment	Temperature range 10 to 40℃ (50 to 104°F) Humidity 20 to 80% RH (Without condensation)
Test Strip	GLUCOCARD X-SENSOR	Dimensions	102x33x18.5mm
Measuring Range	10 TO 600mg/dL (0.6 to 33.3 mmol/L)	Weight	Approx. 30g
Measuring Time	5 seconds	※Design and specifications may be changed without prior notice	
Power source	2 X 3-volt Lithium battery (CR2032)		

arkray global business, inc.

Yousuien-nai, 59 Gansuin-cho,
Kamigyo-ku, Kyoto 602-0008, JAPAN
TEL 81-75-662-8979 FAX +81-75-431-1202