

arkray

Advanced Blood Glucose Test Meter

GLUCOCARD 01-mini plus

GT-1970



More advanced SMBG meter
-User Friendly
-Stylish Design

Advanced Blood Glucose Test Meter

GLUCOCARD 01-mini plus

GT-1970



Design

Impressive & Stylish design

Operation

Comprehensive display
User friendly functions

Small Sample

0.3 μ L micro sampling

Maximum flexibility

Short fill detection
Test strip disposal lever
Test strip slot with light
Backlight LCD



Test strip for GLUCOCARD 01-mini plus

Specification

Test	Glucose in whole blood	Memory	500 test results
Sample Volume	Approx. 0.3 μ L	Operating Environment	Temperature range 10 to 40 $^{\circ}$ C (50 to 104 $^{\circ}$ F) Humidity 20 to 80 % RH (Without condensation)
Test Strip	GLUCOCARD 01-SENSOR	Dimensions	L 102.0 \times W 33.0 \times 18.5 mm
Measuring Range	10 to 600 mg/dL (0.6 to 33.3 mmol/L)	Weight	Approx. 46 g
Measuring Time	7 seconds		
Power source	3 V lithium battery (CR2032) \times 2		

*Design and specifications may be changed without prior notice

ARKRAY GLOBAL BUSINESS, INC.

Yousuien-nai, 59 Gansuin-cho, Kamigyo-ku,
Kyoto 602-0008, JAPAN
TEL +81-75-662-8967 FAX +81-75-662-8973